

— ABLYBIO, Help Your Research



# PPHLN Antibody

货号: **AYP5376**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,THP-1,Mouse testis,Rat testis,Rat spleen,Rat lung
细胞定位	cytosol,Golgi apparatus,nucleoplasm
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human PPHLN.
------	---

## 靶点信息

研究背景	The protein encoded by this gene is one of the several proteins that become sequentially incorporated into the cornified cell envelope during the terminal differentiation of keratinocyte at the outer layers of epidermis. This protein interacts with periplakin, which is known as a precursor of the cornified cell envelope. The cellular localization pattern and insolubility of this protein suggest that it may play a role in epithelial differentiation and contribute to epidermal integrity and barrier formation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]
基因ID	51535
基因名	PPHLN1
Swiss	Q8NEY8 ( <a href="https://www.uniprot.org/uniprotkb/Q8NEY8/entry">https://www.uniprot.org/uniprotkb/Q8NEY8/entry</a> )
别名	CR,HSPC206,HSPC232,PPHLN1,PPHLN Antibody,CDC7 expression repressor,Gastric cancer antigen Ga50,PPHLN

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)