

— ABLYBIO, Help Your Research



# ANKK1 Antibody

货号: **AYP5378**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	cell projection,cytoplasm,nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ANKK1.
------	---

## 靶点信息

研究背景	The protein encoded by this gene belongs to the Ser/Thr protein kinase family, and protein kinase superfamily involved in signal transduction pathways. This gene is closely linked to DRD2 gene (GeneID:1813) on chr 11, and a well studied restriction fragment length polymorphism (RFLP) designated TaqIA, was originally associated with the DRD2 gene, however, later was determined to be located in exon 8 of ANKK1 gene (PMIDs: 18621654, 15146457), where it causes a nonconservative amino acid substitution. It is not clear if this gene plays any role in neuropsychiatric disorders previously associated with TaqIA RFLP.
基因ID	255239
基因名	ANKK1
Swiss	Q8NFD2 ( <a href="https://www.uniprot.org/uniprotkb/Q8NFD2/entry">https://www.uniprot.org/uniprotkb/Q8NFD2/entry</a> )
别名	ANKK1,PKK2,ANKK1 Antibody,Protein kinase PKK2,Sugen kinase 288,X-kinase,SGK288

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))