

— ABLYBIO, Help Your Research



# MUC16/CA125 Antibody

货号: **AYP5401**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF ELISA
推荐浓度	<b>IF:</b> 1:50 - 1:200
理论分子量	1519kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3
细胞定位	extracellular exosome,Golgi lumen,plasma membrane
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MUC16/CA125.
------	---

## 靶点信息

研究背景	This gene encodes a protein that is a member of the mucin family. Mucins are high molecular weight, O-glycosylated proteins that play an important role in forming a protective mucous barrier, and are found on the apical surfaces of the epithelia. The encoded protein is a membrane-tethered mucin that contains an extracellular domain at its amino terminus, a large tandem repeat domain, and a transmembrane domain with a short cytoplasmic domain. The amino terminus is highly glycosylated, while the repeat region contains 156 amino acid repeats unit that are rich in serines, threonines, and prolines. Interspersed within the repeats are Sea urchin sperm protein Enterokinase and Agrin (SEA) modules, leucine-rich repeats and ankyrin (ANK) repeats. These regions together form the ectodomain, and there is a potential cleavage site found near an SEA module close to the transmembrane domain. This protein is thought to play a role in forming a barrier, protecting epithelial cells from pathogens. Products of this gene have been used as a marker for different cancers, with higher expression levels associated with poorer outcomes.
基因ID	94025
基因名	MUC16
Swiss	Q8WXI7 ( <a href="https://www.uniprot.org/uniprotkb/Q8WXI7/entry">https://www.uniprot.org/uniprotkb/Q8WXI7/entry</a> )
别名	CA125,MUC16/CA125 Antibody,MUC16,Ovarian cancer-related tumor marker CA125,Ovarian carcinoma antigen CA125,MUC16/CA125

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)