

# TAF15 Antibody

货号: **AYP5410**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	<a href="#">WB</a> <a href="#">IHC</a> <a href="#">ELISA</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Jurkat
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human TAF15.
------	---

## 靶点信息

研究背景	This gene encodes a member of the TET family of RNA-binding proteins. The encoded protein plays a role in RNA polymerase II gene transcription as a component of a distinct subset of multi-subunit transcription initiation factor TFIID complexes. Translocations involving this gene play a role in acute leukemia and ext raskelatal myxoid chondrosarcoma, and mutations in this gene may play a role in amyotrophic lateral scl erosis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this ge ne.
基因ID	8148
基因名	TAF15

<b>Swiss</b>	Q92804
别名	TAF15;Npl3;RBP56;TAF2N;TAFII68

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)