

— ABLYBIO, Help Your Research



TAF15 Antibody

货号: **AYP5410**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	61kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Jurkat
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TAF15.
------	---

靶点信息

研究背景	This gene encodes a member of the TET family of RNA-binding proteins. The encoded protein plays a role in RNA polymerase II gene transcription as a component of a distinct subset of multi-subunit transcription initiation factor TFIID complexes. Translocations involving this gene play a role in acute leukemia and extr askeletal myxoid chondrosarcoma, and mutations in this gene may play a role in amyotrophic lateral scler osis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gen e.
基因ID	8148
基因名	TAF15
Swiss	Q92804 (https://www.uniprot.org/uniprotkb/Q92804/entry)
别名	TAF15,Npl3,RBP56,TAF2N,TAFI168,TAF15 Antibody,68 kDa TATA-binding protein-associated factor,RNA-bin ding protein 56

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)