

— ABLYBIO, Help Your Research



TP53INP1 Antibody

货号: **AYP5427**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa/27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen,Mouse lung
细胞定位	Cytoplasm,Cytoplasmic vesicle,Nucleus,PML body,autophagosome,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TP53INP1.
------	--

靶点信息

研究背景	Antiproliferative and proapoptotic protein involved in cell stress response which acts as a dual regulator of transcription and autophagy. Acts as a positive regulator of autophagy. In response to cellular stress or activation of autophagy, relocates to autophagosomes where it interacts with autophagosome-associated proteins GABARAP, GABARAPL1/L2, MAP1LC3A/B/C and regulates autophagy. Acts as an antioxidant and plays a major role in p53/TP53-driven oxidative stress response. Possesses both a p53/TP53-independent intracellular reactive oxygen species (ROS) regulatory function and a p53/TP53-dependent transcription regulatory function. Positively regulates p53/TP53 and p73/TP73 and stimulates their capacity to induce apoptosis and regulate cell cycle. In response to double-strand DNA breaks, promotes p53/TP53 phosphorylation on 'Ser-46' and subsequent apoptosis. Acts as a tumor suppressor by inducing cell death by an autophagy and caspase-dependent mechanism. Can reduce cell migration by regulating the expression of SPARC.
基因ID	94241
基因名	TP53INP1
Swiss	Q96A56 (https://www.uniprot.org/uniprotkb/Q96A56/entry)
别名	TP53INP1,SIP,TP53DINP1,TP53INP1A,TP53INP1B,Teap,p53DINP1,TP53INP1 Antibody,Stress-induced protein,p53-dependent damage-inducible nuclear protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)