

ALPK1 Antibody

货号: **AYP5447**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,SW620
细胞定位	centrosome,cilium,cytosol,spindle pole
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ALPK1.
------	---

靶点信息

研究背景	This gene encodes an alpha kinase. Mice which were homozygous for disrupted copies of this gene exhibited coordination defects (PMID: 21208416). Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	80216
基因名	ALPK1
Swiss	Q96QP1
别名	

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn