

— ABLYBIO, Help Your Research



ATF6B Antibody

货号: **AYP5452**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	76kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse lung,mouse brain
细胞定位	Golgi apparatus,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ATF6B.
------	---

靶点信息

研究背景	The protein encoded by this gene is a transcription factor in the unfolded protein response (UPR) pathway during ER stress. Either as a homodimer or as a heterodimer with ATF6-alpha, the encoded protein binds to the ER stress response element, interacting with nuclear transcription factor Y to activate UPR target genes. The protein is normally found in the membrane of the endoplasmic reticulum; however, under ER stress, the N-terminal cytoplasmic domain is cleaved from the rest of the protein and translocates to the nucleus. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	1388
基因名	ATF6B
Swiss	Q99941 (https://www.uniprot.org/uniprotkb/Q99941/entry)
别名	ATF6B,CREB-RP,CREBL1,G13,ATF6B Antibody,Activating transcription factor 6 beta,Protein G13,cAMP response element-binding protein-related protein,cAMP-responsive element-binding protein-like 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)