

— ABLYBIO, Help Your Research



C1orf57 Antibody

货号: **AYP5458**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,A-431,293T
细胞定位	
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human C1orf57.
------	---

靶点信息

研究背景	The protein encoded by this gene is a non-specific nucleoside triphosphatase that is slow-acting in vitro. This gene is overexpressed in many tumor tissues, and while it is not essential for the cell, overexpression is cytotoxic. However, the cytotoxicity is not related to its triphosphatase activity.
------	--

基因ID	84284
基因名	NTPCR
Swiss	Q9BSD7 (https://www.uniprot.org/uniprotkb/Q9BSD7/entry)
别名	NTPCR,C1orf57,HCR-NTPase,C1orf57 Antibody,Nucleoside triphosphate phosphohydrolase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)