

— ABLYBIO, Help Your Research



TINF2 Antibody

货号: **AYP5459**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	15kDa/39kDa/50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HepG2,Mouse skeletal muscle,Mouse heart
细胞定位	Chromosome,Nucleus,Nucleus matrix,telomere
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TINF2.
------	---

靶点信息

研究背景	This gene encodes one of the proteins of the shelterin, or telosome, complex which protects telomeres by allowing the cell to distinguish between telomeres and regions of DNA damage. The protein encoded by this gene is a critical part of shelterin; it interacts with the three DNA-binding proteins of the shelterin complex, and it is important for assembly of the complex. Mutations in this gene cause dyskeratosis congenita (DKC), an inherited bone marrow failure syndrome.
基因ID	26277
基因名	TINF2
Swiss	Q9BSI4 (https://www.uniprot.org/uniprotkb/Q9BSI4/entry)
别名	TINF2,DKCA3,TIN2,TINF2 Antibody,TRF1-interacting nuclear protein 2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)