

— ABLYBIO, Help Your Research



# ILKAP Antibody

货号: **AYP5478**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,Mouse spleen,Mouse testis
细胞定位	Cytoplasm
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ILKAP.
------	---

## 靶点信息

研究背景	The protein encoded by this gene is a protein serine/threonine phosphatase of the PP2C family. This protein can interact with integrin-linked kinase (ILK/ILK1), a regulator of integrin mediated signaling, and regulate the kinase activity of ILK. Through the interaction with ILK, this protein may selectively affect the signaling process of ILK-mediated glycogen synthase kinase 3 beta (GSK3beta), and thus participate in Wnt signaling pathway.
基因ID	80895
基因名	ILKAP
Swiss	Q9H0C8 ( <a href="https://www.uniprot.org/uniprotkb/Q9H0C8/entry">https://www.uniprot.org/uniprotkb/Q9H0C8/entry</a> )
别名	ILKAP,ILKAP2,ILKAP3,PP2C-DELTA,PPM10,ILKAP Antibody

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))