

EKI2 Antibody

货号: **AYP5504**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa/44kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,DU145,HepG2
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human EK12.
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of choline/ethanolamine kinase family which catalyzes the first step of phosphatidylethanolamine (PtdEtn) biosynthesis via the cytidine diphosphate (CDP) ethanolamine pathway. Alternative splicing results in multiple transcript variants.
基因ID	55224
基因名	ETNK2
Swiss	Q9NVF9

别名	ETNK2;EKI2;HMFT1716
----	---------------------

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn