

— ABLYBIO, Help Your Research



# EKI2 Antibody

货号: **AYP5504**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	40kDa/44kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,DU145,HepG2
细胞定位	cytoplasm,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human EKI2.
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of choline/ethanolamine kinase family which catalyzes the first step of phosphatidylethanolamine (PtdEtn) biosynthesis via the cytidine diphosphate (CDP) ethanolamine pathway. Alternative splicing results in multiple transcript variants.
------	--

基因ID	55224
基因名	ETNK2
Swiss	Q9NVF9 ( <a href="https://www.uniprot.org/uniprotkb/Q9NVF9/entry">https://www.uniprot.org/uniprotkb/Q9NVF9/entry</a> )
别名	ETNK2,EKI2,HMFT1716,EKI2 Antibody,Ethanolamine kinase-like protein

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)