

— ABLYBIO, Help Your Research



# MARK2 Antibody

货号: **AYP5516**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	77-87kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Raji,BxPC-3,BT-474,mouse testis,mouse heart,rat testis
细胞定位	Cell membrane,Cytoplasm,Lateral cell membrane,Peripheral membrane protein,cytoskeleton
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human MARK2.
------	---

## 靶点信息

研究背景	This gene encodes a member of the Par-1 family of serine/threonine protein kinases. The protein is an important regulator of cell polarity in epithelial and neuronal cells, and also controls the stability of microtubules through phosphorylation and inactivation of several microtubule-associated proteins. The protein localizes to cell membranes. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	2011
基因名	MARK2
Swiss	Q7KZI7 ( <a href="https://www.uniprot.org/uniprotkb/Q7KZI7/entry">https://www.uniprot.org/uniprotkb/Q7KZI7/entry</a> )
别名	MARK2,EMK-1,EMK1,PAR-1,Par-1b,Par1b,MARK2 Antibody,ELKL motif kinase 1,MAP/microtubule affinity-regulating kinase 2,PAR1 homolog,PAR1 homolog b

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn//www.ablybio.cn>)