

— ABLYBIO, Help Your Research



# KS6A6 Antibody

货号: **AYP5545**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	83kDa/84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human KS6A6.
------	---

## 靶点信息

研究背景	This gene encodes a member of ribosomal S6 kinase family, serine-threonine protein kinases which are regulated by growth factors. The encoded protein may be distinct from other members of this family, however, as studies suggest it is not growth factor dependent and may not participate in the same signaling pathways.
基因ID	27330
基因名	RPS6KA6
Swiss	Q9UK32 ( <a href="https://www.uniprot.org/uniprotkb/Q9UK32/entry">https://www.uniprot.org/uniprotkb/Q9UK32/entry</a> )
别名	RPS6KA6,PP90RSK4,RSK-4,RSK4,S6K-alpha-6,p90RSK6,KS6A6 Antibody,90 kDa ribosomal protein S6 kinase 6,Ribosomal S6 kinase 4,KS6A6

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))