

— ABLYBIO, Help Your Research



# EMCN Antibody

货号: **AYP5550**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	18kDa/26kDa/27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,NIH/3T3
细胞定位	Cell membrane,Membrane,Secreted,Single-pass type I membrane protein,Single-pass type I membrane p rotein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human EMCN.
------	--

## 靶点信息

研究背景	EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al., 2001 [PubMed 11418125]).[supplied by OMIM, Mar 2008]
基因ID	51705
基因名	EMCN
Swiss	Q9ULC0 ( <a href="https://www.uniprot.org/uniprotkb/Q9ULC0/entry">https://www.uniprot.org/uniprotkb/Q9ULC0/entry</a> )
别名	EMCN,EMCN2,MUC14,endomucin,EMCN Antibody,Endomucin-2,Gastric cancer antigen Ga34,Mucin-14

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)