

ICK Antibody

货号: **AYP5558**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF ELISA
推荐浓度	IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	34kDa/71kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SW620,Mouse liver,Mouse lung,Mouse brain
细胞定位	Cell projection,Cytoplasm,Nucleus,cilium,cilium basal body,cytoskeleton,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ICK.
------	---

靶点信息

研究背景	Eukaryotic protein kinases are enzymes that belong to a very extensive family of proteins which share a conserved catalytic core common with both serine/threonine and tyrosine protein kinases. This gene encodes an intestinal serine/threonine kinase harboring a dual phosphorylation site found in mitogen-activating protein (MAP) kinases. The protein localizes to the intestinal crypt region and is thought to be important in intestinal epithelial cell proliferation and differentiation. Alternative splicing has been observed at this locus and two variants, encoding the same isoform, have been identified.
基因ID	22858
基因名	ICK

Swiss	Q9UPZ9
别名	ICK;ECO;LCK2;MRK

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn