

— ABLYBIO, Help Your Research



ZNF148 Antibody

货号: **AYP5560**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	14kDa/88kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Mouse lung,Rat lung
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ZNF148.
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the Kruppel family of zinc finger DNA binding proteins. The encoded protein activates transcription of the T-cell receptor and intestinal alkaline phosphatase genes but represses transcription of the ornithine decarboxylase, vimentin, gastrin, stomelysin, and enolase genes. Increased expression of this gene results in decreased patient survival rates from colorectal cancer, while mutations in this gene have been associated with global developmental delay, hypoplastic corpus callosum, and dysmorphic facies.
基因ID	7707
基因名	ZNF148
Swiss	Q9UQR1 (https://www.uniprot.org/uniprotkb/Q9UQR1/entry)
别名	ZNF148, BERF-1, BFCOL1, GDACCF, HT-BETA, ZBP-89, ZFP148, pHZ-52, ZNF148 Antibody, Transcription factor ZBP-89, Zinc finger DNA-binding protein 89, ZBP89

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)