

— ABLYBIO, Help Your Research



ATF5 Antibody

货号: **AYP5565**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	30kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	rat lung
细胞定位	centrosome,cytosol,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ATF5.
------	--

靶点信息

研究背景	Transcription factor that either stimulates or represses gene transcription through binding of different DNA regulatory elements such as cAMP response element (CRE (consensus: 5'-GTGACGT[AC][AG]-3', ATF5-specific response element (ARE (consensus: 5'-C[CT]TCT[CT]CCTT[AT]-3' but also the amino acid response element (AARE, present in many viral and cellular promoters. Critically involved, often in a cell type-dependent manner, in cell survival, proliferation, and differentiation. Its transcriptional activity is enhanced by CCND3 and slightly inhibited by CDK4. Important regulator of the cerebral cortex formation, functions in cerebral cortical neuroprogenitor cells to maintain proliferation and to block differentiation into neurons. Must be down-regulated in order for such cells to exit the cycle and differentiate (By similarity. Participates in the pathways by which SHH promotes cerebellar granule neuron progenitor cells proliferation (By similarity. Critical for survival of mature olfactory sensory neurons (OSN, directs expression of OSN-specific genes (By similarity.
基因ID	22809
基因名	ATF5
Swiss	Q9Y2D1 (https://www.uniprot.org/uniprotkb/Q9Y2D1/entry)
别名	ATFX,HMFN0395,ATF5,ATF5 Antibody,Activating transcription factor 5,Transcription factor ATFx

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)