

SUPT16H Antibody

货号: **AYP5578**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	119kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HCT116,Jurkat,C6,Mouse brain
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SUPT16H.
------	---

靶点信息

研究背景	Transcription of protein-coding genes can be reconstituted on naked DNA with only the general transcript ion factors and RNA polymerase II. However, this minimal system cannot transcribe DNA packaged into c hromatin, indicating that accessory factors may facilitate access to DNA. One such factor, FACT (facilitate s chromatin transcription), interacts specifically with histones H2A/H2B to effect nucleosome disassembly and transcription elongation. FACT is composed of an 80 kDa subunit and a 140 kDa subunit; this gene en codes the 140 kDa subunit.
基因ID	11198
基因名	SUPT16H

Swiss	Q9Y5B9
别名	CDC68; SPT16; FACTP140; SPT16/CDC68

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn