

— ABLYBIO, Help Your Research



# SPIN1 Antibody

货号: **AYP5583**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,SW480,293T,Jurkat,HeLa,HepG2,SH-SY5Y,Mouse brain
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human SPIN1.
------	---

## 靶点信息

研究背景	Chromatin reader that specifically recognizes and binds histone H3 both trimethylated at 'Lys-4' and asymmetrically dimethylated at 'Arg-8' (H3K4me3 and H3R8me2a) and acts as an activator of Wnt signaling pathway downstream of PRMT2. In case of cancer, promotes cell cancer proliferation via activation of the Wnt signaling pathway. Overexpression induces metaphase arrest and chromosomal instability. Localizes to active rDNA loci and promotes the expression of rRNA genes. May play a role in cell-cycle regulation during the transition from gamete to embryo. Involved in oocyte meiotic resumption, a process that takes place before ovulation to resume meiosis of oocytes blocked in prophase I: may act by regulating maternal transcripts to control meiotic resumption.
基因ID	10927
基因名	SPIN1
Swiss	Q9Y657 ( <a href="https://www.uniprot.org/uniprotkb/Q9Y657/entry">https://www.uniprot.org/uniprotkb/Q9Y657/entry</a> )
别名	SPIN1,SPIN,TDRD24,SPIN1 Antibody,Ovarian cancer-related protein,Spindlin1,OCR

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)