

# 14-3-3 zeta Antibody

货号: **AYP5594**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	19kDa/27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,MCF7,Mouse kidney,Mouse brain,Mouse lung,Rat brain
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human 14-3-3 zeta.
------	---

## 靶点信息

研究背景	This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.
基因ID	7534
基因名	YWHAZ

<b>Swiss</b>	P63104
别名	YWHAZ;14-3-3-zeta;HEL-S-3;HEL-S-93;HEL4;KCIP-1;YWHAD

产品验证

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)