

— ABLYBIO, Help Your Research



14-3-3 zeta Antibody

货号: **AYP5594**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	19kDa/27kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,MCF7,Mouse kidney,Mouse brain,Mouse lung,Rat brain
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human 14-3-3 zeta.
------	---

靶点信息

研究背景	This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.
基因ID	7534
基因名	YWHAZ
Swiss	P63104 (https://www.uniprot.org/uniprotkb/P63104/entry)
别名	YWHAZ,14-3-3-zeta,HEL-S-3,HEL-S-93,HEL4,KCIP-1,YWHAD,14-3-3 zeta Antibody,Protein kinase C inhibitor protein 1,14-3-3 zeta

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)