

— ABLYBIO, Help Your Research



ADCY1 Antibody

货号: **AYP5606**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	123kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse brain,mouse kidney
细胞定位	Cell membrane,Cytoplasm,Membrane,Membrane raft,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ADCY1.
------	---

靶点信息

研究背景	This gene encodes a member of the of adenylate cyclase gene family that is primarily expressed in the brain. This protein is regulated by calcium/calmodulin concentration and may be involved in brain development. Alternate splicing results in multiple transcript variants.
基因ID	107
基因名	ADCY1
Swiss	Q08828 (https://www.uniprot.org/uniprotkb/Q08828/entry)
别名	ADCY1,AC1,DFNB44,ADCY1 Antibody,ATP pyrophosphate-lyase 1,Adenylate cyclase type I,Adenylyl cyclase 1,Ca(2+)/calmodulin-activated adenylyl cyclase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)