

— ABLYBIO, Help Your Research



IFIT1 Antibody

货号: **AYP5684**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/55kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human IFIT1.
------	-----------------------------------------------

靶点信息

研究背景	This gene encodes a protein containing tetratricopeptide repeats that was originally identified as induced upon treatment with interferon. The encoded protein may inhibit viral replication and translational initiation. This gene is located in a cluster on chromosome 10 with five other closely related genes. There is a pseudogene for this gene on chromosome 13. Alternatively spliced transcript variants encoding multiple isoforms have been observed.
基因ID	3434
基因名	IFIT1
Swiss	P09914 (https://www.uniprot.org/uniprotkb/P09914/entry)
别名	IFIT1,C56,G10P1,IFI-56,IFI-56K,IFI56,IFIT-1,IFNAI1,ISG56,P56,RNM561,IFIT1 Antibody,Interferon-induced 56 kDa protein,Interferon-induced protein with tetratricopeptide repeats 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)