

— ABLYBIO, Help Your Research



Galectin 1 Antibody

货号: **AYP5695**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS,A-549,Mouse testis,Mouse thymus,Rat lung,Rat heart
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human Galectin 1.
------	--

靶点信息

研究背景	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. This gene product may act as an autocrine negative growth factor that regulates cell proliferation.
------	--

基因ID	3956
基因名	LGALS1
Swiss	P09382 (https://www.uniprot.org/uniprotkb/P09382/entry)
别名	LGALS1,GAL1,GBP,Galectin 1 Antibody,14 kDa laminin-binding protein,14 kDa lectin,Beta-galactoside-binding lectin L-14-I,Galaptin,HBL,HPL,Lactose-binding lectin 1,Lectin galactoside-binding soluble 1,Putative MAPK-activating protein PM12,S-Lac lectin 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)