

— ABLYBIO, Help Your Research



MRPL11 Antibody

货号: **AYP5722**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF: 1:50 - 1:200
理论分子量	18kDa/19kDa/20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,HeLa,MCF7,A-549,Mouse kidney,Mouse brain,Mouse heart,Rat heart
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MRPL11.
------	--

靶点信息

研究背景	This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12.
基因ID	65003
基因名	MRPL11
Swiss	Q9Y3B7 (https://www.uniprot.org/uniprotkb/Q9Y3B7/entry)
别名	MRPL11,CGI-113,L11MT,MRP-L11,MRPL11 Antibody,39S ribosomal protein L11,mitochondrial

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)