

— ABLYBIO, Help Your Research



MRPL18 Antibody

货号: **AYP5728**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,K-562
细胞定位	extracellular space,mitochondrial inner membrane,mitochondrial ribosome,mitochondrion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human MRPL18.
------	--

靶点信息

研究背景	This nuclear gene encodes a protein component of the larger 39S subunit of mitochondrial ribosome. This protein may also aid in the import of nuclear-encoded 5S rRNA into mitochondria. Alternative splicing results in multiple transcript variants, most of which are not predicted to encode a protein. A pseudogene of this gene is found on chromosome 16.
基因ID	29074
基因名	MRPL18
Swiss	Q9H0U6 (https://www.uniprot.org/uniprotkb/Q9H0U6/entry)
别名	MRPL18,HSPC071,L18mt,MRP-L18,MRPL18 Antibody,39S ribosomal protein L18,mitochondrial

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)