

— ABLYBIO, Help Your Research



RPL14 Antibody

货号: **AYP5762**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	23kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,HeLa,293T,Jurkat,HepG2,Mouse ovary
细胞定位	cytoplasm,cytosol,cytosolic ribosome,extracellular exosome,postsynaptic density
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human RPL14.
------	---

靶点信息

研究背景	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L14E family of ribosomal proteins. It contains a basic region-leucine zipper (bZIP)-like domain. The protein is located in the cytoplasm. This gene contains a trinucleotide (GCT) repeat tract whose length is highly polymorphic; these triplet repeats result in a stretch of alanine residues in the encoded protein. Transcript variants utilizing alternative polyA signals and alternative 5'-terminal exons exist but all encode the same protein. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.
基因ID	9045
基因名	RPL14
Swiss	P50914 (https://www.uniprot.org/uniprotkb/P50914/entry)
别名	RPL14,CAG-ISL-7,CTG-B33,L14,RL14,hRL14,RPL14 Antibody,60S ribosomal protein L14,CAG-ISL 7

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)