

— ABLYBIO, Help Your Research



# RPL17 Antibody

货号: **AYP5764**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	17kDa/21kDa/26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293T,A-549,SW480,HeLa,Mouse lung,Mouse thymus,Rat liver
细胞定位	cytoplasm,cytosol,cytosolic ribosome,nucleus
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human RPL17.
------	---

## 靶点信息

研究背景	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L22P family of ribosomal proteins. It is located in the cytoplasm. This gene has been referred to as rpL23 because the encoded protein shares amino acid identity with ribosomal protein L23 from <i>Halobacterium marismortui</i> ; however, its official symbol is RPL17. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream C18orf32 (chromosome 18 open reading frame 32) gene.
基因ID	6139
基因名	RPL17
Swiss	P18621 ( <a href="https://www.uniprot.org/uniprotkb/P18621/entry">https://www.uniprot.org/uniprotkb/P18621/entry</a> )
别名	RPL17,L17,PD-1,RPL23,RPL17 Antibody,60S ribosomal protein L17,60S ribosomal protein L23

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))