

— ABLYBIO, Help Your Research



# ACTR3 Antibody

货号: **AYP5794**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF:</b> 1:50 - 1:200
理论分子量	47kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,HL-60,HeLa,Jurkat,COS-1,COS-7,Mouse spleen
细胞定位	Cell projection,Cytoplasm,cytoskeleton
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ACTR3.
------	---

## 靶点信息

研究背景	The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Three transcript variants encoding two different isoforms have been found for this gene.
基因ID	10096
基因名	ACTR3
Swiss	P61158 ( <a href="https://www.uniprot.org/uniprotkb/P61158/entry">https://www.uniprot.org/uniprotkb/P61158/entry</a> )
别名	ACTR3,ARP3,ACTR3 Antibody,Actin-like protein 3

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)