

— ABLYBIO, Help Your Research



ACOT2 Antibody

货号: **AYP5799**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa/53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,U-87MG,Mouse kidney
细胞定位	Mitochondrion
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ACOT2.
------	-----------------------------------------------

靶点信息

研究背景	This gene encodes a member of the acyl-CoA thioesterase protein family, and is one of four acyl-CoA hydrolase genes located in a cluster on chromosome 14. Alternative splicing results in multiple transcript variants.
------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

基因ID	10965
基因名	ACOT2
Swiss	P49753 (https://www.uniprot.org/uniprotkb/P49753/entry)
别名	ACOT2,CTE-IA,CTE1A,MTE1,PTE2,PTE2A,ZAP128,ACOT2 Antibody,Acyl-coenzyme A thioester hydrolase 2 a,Long-chain acyl-CoA thioesterase 2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)