

— ABLYBIO, Help Your Research



SLC25A6 Antibody

货号: **AYP5806**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	32kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,293T,K-562,Mouse thymus
细胞定位	Mitochondrion inner membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human SLC25A6.
------	---

靶点信息

研究背景	This gene is a member of the mitochondrial carrier subfamily of solute carrier protein genes. The product of this gene functions as a gated pore that translocates ADP from the cytoplasm into the mitochondrial matrix and ATP from the mitochondrial matrix into the cytoplasm. The protein is implicated in the function of the permeability transition pore complex (PTPC), which regulates the release of mitochondrial products that induce apoptosis. The human genome contains several non-transcribed pseudogenes of this gene.
基因ID	293
基因名	SLC25A6
Swiss	P12236 (https://www.uniprot.org/uniprotkb/P12236/entry)
别名	SLC25A6,AAC3,ANT,ANT 2,ANT 3,ANT3,ANT3Y,ANT2,SLC25A6 Antibody,ADP,ATP carrier protein 3,ATP carrier protein,isoform T2,Adenine nucleotide translocator 3,Solute carrier family 25 member 6

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)