

— ABLYBIO, Help Your Research



ARF4 Antibody

货号: **AYP5809**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse brain,Mouse testis,Rat testis,Rat brain
细胞定位	Golgi apparatus,Lipid-anchor,Membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human ARF4.
------	--

靶点信息

研究背景	<p>This gene is a member of the human ARF gene family whose members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The gene products include 5 ARF proteins and 11 ARF-like proteins and constitute one family of the RAS superfamily. The ARF proteins are categorized as class I, class II and class III; this gene is a class II member. The members of each class share a common gene organization. The ARF4 gene spans approximately 12kb and contains six exons and five introns. This gene is the most divergent member of the human ARFs. Conflicting map positions at 3p14 or 3p21 have been reported for this gene.</p>
基因ID	378
基因名	ARF4
Swiss	P18085 (https://www.uniprot.org/uniprotkb/P18085/entry)
别名	ARF4,ARF2,ARF4 Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)