

— ABLYBIO, Help Your Research



AIBP Antibody

货号: **AYP5845**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	20kDa/31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,NIH/3T3,Rat lung
细胞定位	Mitochondrion,Secreted
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human AIBP.
------	--

靶点信息

研究背景	The product of this gene interacts with apolipoprotein A-I (apoA-I), the major apolipoprotein of high-density lipoproteins (HDLs). It is secreted into some bodily fluids, and its synthesis and secretion are stimulated in vitro by incubating cells with apoA-I. The human genome contains related pseudogenes.
基因ID	128240
基因名	NAXE
Swiss	Q8NCW5 (https://www.uniprot.org/uniprotkb/Q8NCW5/entry)
别名	NAXE,AIBP,APOA1BP,PEBEL,YJEFN1,AIBP Antibody,Apolipoprotein A-I-binding protein,NAD(P)HX epimerase,YjeF N-terminal domain-containing protein 1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)