

— ABLYBIO, Help Your Research



# ARSA Antibody

货号: **AYP5856**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF:</b> 1:50 - 1:200
理论分子量	44kDa/53kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,HepG2,A-549,Mouse testis,Mouse liver,Mouse brain
细胞定位	Lysosome
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human ARSA.
------	--

## 靶点信息

研究背景	The protein encoded by this gene hydrolyzes cerebroside sulfate to cerebroside and sulfate. Defects in this gene lead to metachromatic leucodystrophy (MLD), a progressive demyelination disease which results in a variety of neurological symptoms and ultimately death. Alternatively spliced transcript variants have been described for this gene.
基因ID	410
基因名	ARSA
Swiss	P15289 ( <a href="https://www.uniprot.org/uniprotkb/P15289/entry">https://www.uniprot.org/uniprotkb/P15289/entry</a> )
别名	ARSA,MLD,ARSA Antibody,Cerebroside-sulfatase

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)