

— ABLYBIO, Help Your Research



DDX51 Antibody

货号: **AYP5883**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	72kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Jurkat,Mouse brain,Mouse spleen,Rat testis
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human DDX51.
------	---

靶点信息

研究背景	ATP-binding RNA helicase involved in the biogenesis of 60S ribosomal subunits.
基因ID	317781

基因名	DDX51
Swiss	Q8N8A6 (https://www.uniprot.org/uniprotkb/Q8N8A6/entry)
别名	DDX51,DDX51 Antibody,DEAD box protein 51

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)