

— ABLYBIO, Help Your Research



BCLAF1 Antibody

货号: **AYP5893**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	85kDa/100kDa/105kDa/106kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y
细胞定位	Cytoplasm,Nucleus,Nucleus speckle,nucleoplasm
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human BCLAF1.
------	--

靶点信息

研究背景	This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	9774
基因名	BCLAF1
Swiss	Q9NYF8 (https://www.uniprot.org/uniprotkb/Q9NYF8/entry)
别名	BCLAF1,BTF,bK211L9.1,BCLAF1 Antibody,BCLAF1 and THRAP3 family member 1,KIAA0164

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)