

— ABLYBIO, Help Your Research



B4GALT5 Antibody

货号: **AYP5906**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Golgi apparatus
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human B4GALT5.
------	---

靶点信息

研究背景	This gene is one of seven beta-1,4-galactosyltransferase (beta4GalT) genes. They encode type II membrane-bound glycoproteins that appear to have exclusive specificity for the donor substrate UDP-galactose; a II transfer galactose in a beta1,4 linkage to similar acceptor sugars: GlcNAc, Glc, and Xyl. Each beta4GalT has a distinct function in the biosynthesis of different glycoconjugates and saccharide structures. As type II membrane proteins, they have an N-terminal hydrophobic signal sequence that directs the protein to the Golgi apparatus and which then remains uncleaved to function as a transmembrane anchor. By sequence similarity, the beta4GalTs form four groups: beta4GalT1 and beta4GalT2, beta4GalT3 and beta4GalT4, beta4GalT5 and beta4GalT6, and beta4GalT7. The function of the enzyme encoded by this gene is not clear. This gene was previously designated as B4GALT4 but was renamed to B4GALT5. In the literature it is also referred to as beta4GalT2. [provided by RefSeq, Jul 2008]
基因ID	9334
基因名	B4GALT5
Swiss	O43286 (https://www.uniprot.org/uniprotkb/O43286/entry)
别名	B4Gal-T5,BETA4-GALT-IV,beta4Gal-T5,beta4GalT-V,gt-V,B4GALT5,B4GALT5 Antibody,Beta-1,4-GalT II,Gluco sylceramide beta-1,4-galactosyltransferase,Lactosylceramide synthase,UDP-Gal:beta-GlcNAc beta-1,4-gal actosyltransferase 5,UDP-galactose:beta-N-acetylglucosamine beta-1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)