

— ABLYBIO, Help Your Research



BMP10 Antibody

货号: **AYP5911**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	rat heart
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human BMP10.
------	---

靶点信息

研究背景	The protein encoded by this gene is a member of the TGF-beta family of growth factors. Data suggest that the similar protein in mouse plays an important role in trabeculation of the embryonic heart. In human, this protein may signal through receptor serine/threonine kinases.
------	---

基因ID	27302
基因名	BMP10
Swiss	O95393 (https://www.uniprot.org/uniprotkb/O95393/entry)
别名	BMP10,BMP10 Antibody

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)