

— ABLYBIO, Help Your Research



CHP Antibody

货号: **AYP5928**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	22kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,293T,HeLa,Mouse brain,Mouse heart,Mouse kidney,Rat brain
细胞定位	Cell membrane,Cytoplasm,Endoplasmic reticulum,Endoplasmic reticulum-Golgi intermediate compartment,Lipid-anchor,Membrane,Nucleus,cytoskeleton
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human CHP.
------	---

靶点信息

研究背景	This gene encodes a phosphoprotein that binds to the Na ⁺ /H ⁺ exchanger NHE1. This protein serves as a n essential cofactor which supports the physiological activity of NHE family members and may play a role in the mitogenic regulation of NHE1. The protein shares similarity with calcineurin B and calmodulin and i t is also known to be an endogenous inhibitor of calcineurin activity.
基因ID	11261
基因名	CHP1
Swiss	Q99653 (https://www.uniprot.org/uniprotkb/Q99653/entry)
别名	CHP1,CHP,SLC9A1BP,Sid470p,p22,p24,CHP Antibody,Calcineurin B-like protein,Calcium-binding protein C HP,Calcium-binding protein p22,EF-hand calcium-binding domain-containing protein p22

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)