

— ABLYBIO, Help Your Research



CEP70 Antibody

货号: **AYP5961**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	25kDa/64kDa/67kDa/69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,A-549,Rat testis
细胞定位	Cytoplasm,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human CEP70.
------	---

靶点信息

研究背景	Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.
------	--

基因ID	80321
基因名	CEP70
Swiss	Q8NHQ1 (https://www.uniprot.org/uniprotkb/Q8NHQ1/entry)
别名	CEP70,BITE,CEP70 Antibody,p10-binding protein

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)