

— ABLYBIO, Help Your Research



# CKLF3 Antibody

货号: **AYP5981**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB ELISA
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	6kDa/19kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	231
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息

抗原信息	Synthesized peptide derived from Human CKLF3.
------	---

## 靶点信息

研究背景	This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and the transmembrane 4 superfamilies of signaling molecules. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 16. Alternatively spliced transcript variants containing different 5' UTRs, but encoding the same protein, have been identified.
基因ID	123920
基因名	CMTM3
Swiss	Q96MX0 ( <a href="https://www.uniprot.org/uniprotkb/Q96MX0/entry">https://www.uniprot.org/uniprotkb/Q96MX0/entry</a> )
别名	CMTM3,BNAS2,CKLFSF3,CKLF3 Antibody,Chemokine-like factor superfamily member 3,CKLF3

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)