

— ABLYBIO, Help Your Research



COPZ1 Antibody

货号: **AYP5984**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	14kDa/17kDa/18kDa/20kDa/21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,LO2,HeLa,Mouse kidney,Mouse testis,Mouse spleen
细胞定位	COPI-coated vesicle membrane,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Golgi apparatus membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human COPZ1.
------	---

靶点信息

研究背景	This gene encodes a subunit of the cytoplasmic coatamer protein complex, which is involved in autophagy and intracellular protein trafficking. The coatamer protein complex is comprised of seven subunits and functions as the coat protein of coat protein complex (COP)I-vesicles. Alternative splicing results in multiple transcript variants.
基因ID	22818
基因名	COPZ1
Swiss	P61923 (https://www.uniprot.org/uniprotkb/P61923/entry)
别名	COPZ1,CGI-120,COPZ,HSPC181,zeta-COP,zeta1-COP,COPZ1 Antibody,Zeta-1-coat protein

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)