

— ABLYBIO, Help Your Research



C1S Antibody

货号: **AYP5992**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC ELISA
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	76kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	blood microparticle,complement component C1 complex,extracellular region,extracellular space
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human C1S.
------	---

靶点信息

研究背景	This gene encodes a serine protease, which is a major constituent of the human complement subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are the cause of selective C1s deficiency.
基因ID	716
基因名	C1S
Swiss	P09871 (https://www.uniprot.org/uniprotkb/P09871/entry)
别名	C1S,EDSPD2,C1S Antibody,C1 esterase,Complement component 1 subcomponent s,Complement C1s subcomponent chain A,Complement C1s subcomponent chain B

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)