

— ABLYBIO, Help Your Research



DBH Antibody

货号: **AYP6038**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	69kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cytoplasmic vesicle,Single-pass type II membrane protein,chromaffin granule lumen,chromaffin granule membrane,secretory vesicle,secretory vesicle lumen,secretory vesicle membrane
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human DBH.
------	---

靶点信息

研究背景	The protein encoded by this gene is an oxidoreductase belonging to the copper type II, ascorbate-dependent monooxygenase family. It is present in the synaptic vesicles of postganglionic sympathetic neurons and converts dopamine to norepinephrine. It exists in both soluble and membrane-bound forms, depending on the absence or presence, respectively, of a signal peptide.
基因ID	1621
基因名	DBH
Swiss	P09172 (https://www.uniprot.org/uniprotkb/P09172/entry)
别名	DBH,DBM,DBH Antibody,Dopamine beta-monoxygenase

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)