

— ABLYBIO, Help Your Research



CCNB1IP1 Antibody

货号: **AYP6041**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	31kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	mouse heart
细胞定位	Nucleus. Chromosome
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human CCNB1IP1.
------	--

靶点信息

研究背景	HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.
------	--

基因ID	57820
基因名	CCNB1IP1
Swiss	Q9NPC3 (https://www.uniprot.org/uniprotkb/Q9NPC3/entry)
别名	CCNB1IP1,C14orf18,HEI10,CCNB1IP1 Antibody,Cyclin-B1-interacting protein 1,Human enhancer of invasion 10,RING-type E3 ubiquitin transferase CCNB1IP1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)