

— ABLYBIO, Help Your Research



TUBGCP3 Antibody

货号: **AYP6063**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC ELISA
推荐浓度	IHC: 1:50 - 1:200
理论分子量	48kDa/93kDa/103kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse spleen,Mouse thymus,Rat testis
细胞定位	Cytoplasm,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TUBGCP3.
------	---

靶点信息

研究背景	Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.
基因ID	10426

基因名	TUBGCP3
Swiss	Q96CW5 (https://www.uniprot.org/uniprotkb/Q96CW5/entry)
别名	TUBGCP3,104p,GCP3,Grip104,SPBC98,Spc98,Spc98p,TUBGCP3 Antibody,Gamma-ring complex protein 104 kDa,Spindle pole body protein Spc98 homolog

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)