

— ABLYBIO, Help Your Research



TUBGCP5 Antibody

货号: **AYP6064**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF ELISA
推荐浓度	WB: 1:500 - 1:2000 IF: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat heart
细胞定位	centrosome,cytosol,equatorial microtubule organizing center,gamma-tubulin complex,spindle pole
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human TUBGCP5.
------	---

靶点信息

研究背景	Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.
------	--

基因ID	114791
基因名	TUBGCP5
Swiss	Q96RT8 (https://www.uniprot.org/uniprotkb/Q96RT8/entry)
别名	GCP5,TUBGCP5,TUBGCP5 Antibody,KIAA1899

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)