

— ABLYBIO, Help Your Research



USP53 Antibody

货号: **AYP6095**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC ELISA
推荐浓度	IHC: 1:50 - 1:200
理论分子量	120kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,293T,A-431
细胞定位	bicellular tight junction,cell-cell junction
纯化	Affinity purification

抗原信息

抗原信息	Synthesized peptide derived from Human USP53.
------	---

靶点信息

研究背景	Tight junction-associated protein that is involved in the survival of auditory hair cells and hearing. Maybe by modulating the barrier properties and mechanical stability of tight junctions (By similarity. Has no peptidase activity.
------	--

基因ID	54532
基因名	USP53
Swiss	Q70EK8 (https://www.uniprot.org/uniprotkb/Q70EK8/entry)
别名	USP53, USP53 Antibody, Ubiquitin-specific peptidase 53, KIAA1350

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)